

Data Sheet

Im7 Resin (50% suspension)

SKU No.: 01-1701, 01-1702, 01-1705, 01-1715, 01-2701, 01-2705, 01-3701, 01-375

Description	Binds the CL7 tag (~16 kDA) fused to target proteins. Following cleavage with the appropriate protease, the target protein releases from the Im7-bound CL7.
Particle Size	Crosslinked agarose 6B beads (45 - 165 μ M)
pH Stability	The Im7 protein is stable on the beads at pH 3-10 (normal working conditions). However, CL7/Im7 binding is stable at pH 4.2 - 10 only.
Salt stability	\leq 4 M NaCl tested
Binding Capacity	35-40 mg CL7/mL resin
Storage/Shipping Concentration	50/50 buffer/resin slurry. 1-, 2-, 5-, and 15-mL of settled resin in 2, 4, 10, and 30 mL of slurry, respectively.
Recommended Operating Temperature	4° C or room temperature
Reactivation Details	To remove the CL7 protein and reactivate the resin, wash the column with guanidine hydrochloride, exchanging it into physiological buffer.

Additional Information: Target protein characteristics (e.g. protein size, conformation, and concentration); flow rate (i.e. lower flow rates may increase the binding capacity); and other parameters (e.g. pH and temperature) can affect the column binding capacity.

You can download the full protocol from <https://trialtusbioscience.com/products/#protocols>.

**For research use only.
Not for diagnostic or therapeutic use.**

Licensing Information

TriAltus Bioscience holds the exclusive, worldwide license to the CL7 protein purification technology platform. It was licensed from the University of Alabama at Birmingham (UAB) in Birmingham, Alabama, USA. An international patent filing has been made with protection being sought in the United States, Europe, and other major markets. The CL7 purification technology is available for research use. For commercial use or resale, contact us at sales@trialtusbioscience.com to discuss commercial licensing.

Trademarks

The Company name, the terms "TriAltus" and "TriAltus Bioscience", the Company logo, and all related names, logos, product and service names, designs, and slogans are trademarks of the Company or its affiliates or licensors. You must not use such marks without the prior written permission of the Company. All other names, logos, product and service names, designs, and slogans on this Website are the trademarks of their respective owners.

Version: 1.0
Revision Date: 11/20/2019